PRODIGY VOICE[®]

Blood Glucose Monitoring System

Owner's Manual

Dear Prodigy® Meter Owner:

Thank you for choosing the Prodigy Voice® meter as your blood glucose monitoring system. Please read this manual carefully as it contains important information about your new Prodigy® system. A warranty registration card is included with your system. **Please return the completed card to us or complete warranty online at www.prodigymeter.com.**

Prodigy® meters are designed to help you and your healthcare professionals monitor your blood glucose levels. This owner's manual will help you learn how to use your Prodigy® meter effectively. Customer Care is available 24 hours, 7 days a week by calling toll free **1.800.243.2636**. All questions concerning test results should be directed to your healthcare professional.

Prodigy® meters have the latest technology for blood glucose monitoring and are made with quality components. All Prodigy® meters are easy to use, give you fast, accurate test results with a minimal sample of blood, have large, easy-to-read display screens, and are lightweight and portable for your convenience.

Prodigy® meters:

- **Require No Coding**, allowing you to save time and avoid human error due to improper coding.
- Allow you to perform Alternate Site Testing (AST).
- Have memory and data management capabilities. Prodigy's® **free software** gives you and your healthcare professionals powerful graphic tools to manage your diabetes.

Table of Contents

Letter to Meter Owner.....	2-3
Table of Contents	4-5
Important Safety Instructions	6-8
Important Health-Related Information	9-10
Alternate Site Testing (AST).....	11-13
About Prodigy® Test Strips	14-16
Key Functions of the Test Strip.....	14
Important Test Strip Information.....	15-16
Introduction to Prodigy® Meters	17-18
Intended Use.....	17
Test Principle	17-18
Contents of the Prodigy Voice® Meter	18
Setting up the Prodigy Voice® Meter.....	19-30
Key Functions of the Meter.....	19-20
Meter Display	21-22
Speaking Function.....	23-24
Setup Steps.....	25-30
Control Solution Testing.....	31-38
Performing a Control Solution Test	34-38
Performing a Blood Glucose Test	39-51
Preparing the Lancing Device.....	39-41
Getting a Blood Sample	41-43
Tips for Applying Blood Without Sight.....	44-47
Applying Blood to the Test Strip.....	48-51

Using the Meter Memory.....	52-54
Specifications.....	55
Viewing Results on a Personal Computer	56
Caring for your Prodigy® Meter.....	57-61
Cleaning your System	57
Lancing Device	57
Disinfecting your System	57-58
Storing your System	58
Battery	59
Low Battery.....	59-60
Battery Replacement.....	60-61
System Troubleshooting	62-68
Special Messages	62-63
Error Messages	63-64
Troubleshooting Guide.....	65-68
Information About your Prodigy® Meter.....	69-77
Comparing Meter and Laboratory Results.....	69-71
Performance Characteristics	71
Expected Test Results.....	72
Unexpected Test Results.....	73-74
Accuracy of Meter Measurements	75-77
Symbols Information	78
Warranty Information	79
Limited Lifetime Warranty	79

Important Safety Instructions

Read this before using your Prodigy® meter. The following basic safety precautions should always be taken.

- Your Prodigy® meter, test strips, lancets and lancing device are for single patient use only. Do not share them with anyone, including other family members. Do not use on multiple patients.
- Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
- Use the device only for the intended use described in this manual.
- Do not use test strips and control solutions with your Prodigy® meter that are not supplied by Prodigy®.
- Do not use the device if it is not working properly, or if it has suffered any damage.
- Before using any product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

- Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard.
- Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.
- All parts of this kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.
- These devices are intended to be used for patient self-monitoring and should not be used to collect blood from more than one person as this poses a risk of transmitting blood-borne pathogens such as Hepatitis B or HIV.

For further information, please refer to:

"FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Blood-borne Pathogens: Initial Communication" (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

“CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Blood-borne Pathogens” (2012)
<http://cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

DO NOT CHANGE YOUR TREATMENT BASED ON A SINGLE RESULT THAT DOES NOT MATCH HOW YOU FEEL OR IF YOU BELIEVE THAT YOUR TEST RESULT COULD BE INCORRECT.

Important Health-Related Information

Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.

Under normal blood concentrations and conditions, the following should not significantly affect results:

- Elevated blood triglycerides
- Reducing substances such as uric acid and ascorbic acid
- Acetaminophen
- Dopa
- Methyl dopa
- L-dopa
- Tolbutamide

If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, contact your healthcare professional.

Test results below 60 mg/dL (3.3 mmol/L) indicate low blood glucose (hypoglycemia). Test results greater than 240 mg/dL (13.3 mmol/L) indicate high blood glucose (hyperglycemia). If your results are below 60 mg/dL or above 240 mg/dL, repeat the test, and if the results are still below 60 mg/dL (3.3 mmol/L) or above 240 mg/dL (13.3 mmol/L), consult your healthcare professional immediately.

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state. Please refer to the test strip package insert for additional important information.

Alternate Site Testing (AST)

There are important limitations to Alternate Site Testing (AST). Please consult your healthcare professional before performing AST.

What is AST?

Alternate Site Testing (AST) means you can use parts of the body other than your fingertips to check your blood glucose levels. Prodigy® meters allow you to test on your palm, forearm, upper arm, calf or thigh. See Figure 1.


Figure 1

What is the Advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, nerve endings are not so numerous and you will not feel as much pain as you will experience at the fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Therefore, if you are testing your blood glucose level during or immediately after a meal, physical exercise or stressful event, take the blood sample from your fingertip only.

Use AST only:

- ❖ Two hours or more after your last meal.
- ❖ Two hours or more after taking insulin.
- ❖ Two hours or more after exercise.
- ❖ During a relaxed and calm state.

Do not use AST if:

- ❖ You have reason to believe you have hypoglycemia or hyperglycemia.
- ❖ Your routine glucose results are often fluctuating.
- ❖ You are pregnant.

**To increase accuracy when using AST, rub the puncture site before extracting blood. An alternate site device is required to perform AST. If you do not have an alternate site device but would like to use one, please contact Customer Care at 1.800.243.2636.*

About Prodigy® Test Strips

Key Functions of the Test Strip

All Prodigy® meters measure the amount of sugar in the blood, commonly referred to as blood glucose. Blood is applied to the opening of the absorbent channel on the end of the test strip and is automatically drawn into the test strip.

The test strip consists of the following parts:

Absorbent Channel

Apply a drop of blood and it will be drawn in automatically.

Confirmation Window

Shows whether enough blood has been drawn into the test strip's absorbent channel.

Test Strip Handle

Hold this part to insert the test strip into the test strip port on the meter.

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will not go any further.

Figure 2

Please refer to the "Performing a Blood Glucose Test" section for complete instructions.

Important Test Strip Information

- Store test strip packages in a cool, dry place between 39°F- 104°F (4°C-40°C). Keep away from direct sunlight and heat. Do not refrigerate. Do not store in your car.
- Store your test strips in their original vial only. Do not transfer them to a new vial or any other container.
- With clean, dry hands you may touch the test strip anywhere on its surface when removing it from the vial or inserting it into the meter.
- Do not clean the testing site with alcohol. Instead, wash the test site thoroughly with soap and water only.
- Immediately use a test strip after removing it from the vial, replace the vial cap and close it tightly to prevent exposure to moisture.
- Only apply a blood sample or a control

solution sample to the test strip's absorbent channel. Applying other substances to the test strip's absorbent channel will cause inaccurate results.

- Record the discard date on the vial label when you first open it. **Discard remaining test strips 90 days after the first opening date.**
- Do not use test strips beyond the expiration date printed on the package.
- After removing the test strip from the vial, insert it directly into the meter. Do not lay the test strip on any other surface.

Warning: Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard. Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.

Introduction to Prodigy® Meters

Intended Use

Prodigy® blood glucose monitoring systems are intended for use outside the body (in vitro diagnostic use only) and should only be used for blood glucose (blood sugar) testing with fresh capillary whole blood samples.

These systems are intended for use in the home and in clinical settings and should not be used for the diagnosis of diabetes or for the testing of newborns. Consult your healthcare provider if unusual readings occur.

Test Principle

Blood glucose is measured by an electrical current that is produced when a blood sample mixes with the reagent (special chemicals) of the test strip. The electrical current changes with the amount of glucose in the blood sample. Prodigy® meters measure the strength of the electrical current, calculate your blood glucose level and then display your result in either mg/dL or mmol/L.

PRODIGY Voice® Blood Glucose Meter

Prodigy® meters, test strips and control solutions have been designed, tested and proven to work together as a system to produce accurate blood glucose test results.

Important: Use only Prodigy® control solutions and test strips with your Prodigy® meter. Using other test strips and control solutions with this meter can produce inaccurate results.

Contents of the Prodigy Voice® Meter

The Prodigy Voice® meter is available as a meter kit only. Please review the contents of your purchase to confirm that all the components are included as listed on the side of your meter box.

PRODIGY Voice® Blood Glucose Meter

Setting up the Prodigy Voice® Meter

Key Functions of the Meter

Test Strip Slot

Insert the test strip here. The meter will turn on automatically.

LCD Display

Guides you through the test using symbols and simple messages.

Repeat Button

Repeats last message spoken.

Set Button

Used to setup the meter.

Data Port

Port for USB cable connection to your computer.



Figure 3

Key Functions of the Meter

Meter Display

Blood Drop Symbol
Flashes when sample should be applied.

Test Strip Symbol
Appears when the meter is in testing mode.

Test Result Area
Displays test results.


Day Average
Indicates that the display test result is an average.

Date

Time

Control Solution Test Symbol

Shows that you are in control solution mode.


Figure 6

Temperature Symbol

Appears when temperature is outside of operating range.

Voice Symbol

Shows if audio function is on.

Unit of Measure

Indicates what unit of measure the test result is displayed in.

Low Battery Symbol

Appears when the battery power is low.

Memory Symbol

Appears when you review the memory.

Speaking Function

The Prodigy Voice® meter talks you through each step of your setup, testing procedure and memory recall.

When does the meter speak?

When the meter is turned on.

What does the meter say?

Thank you for using Prodigy Voice®. Today's date November 18, 2013. The time 10:08 PM.

When the meter is ready to test (a drop of blood appears on display).

Please apply blood to the test strip

When the test is completed (the result appears on the display).

Your blood glucose (number) milligrams per deciliter/millimoles per liter.

Speaking Function

When does the meter speak?	What does the meter say?
When you turn off the meter	Goodbye.
When a used test strip is inserted	Test strip has been used. Replace with a new one.
When blood is applied to the test strip	Now testing.

Setup Steps

The Prodigy Voice® meter has several different settings that may need to be adjusted before using the meter. You should check the time and date each time the batteries are replaced.

STEP 1: Insert the Batteries

Open the battery cover located on the backside of the meter. Insert two (2) AAA Alkaline batteries and align them correctly. When batteries are installed correctly, the meter will turn on. See Figure 7.


Figure 7

Voice speaks: **"Your Prodigy meter is on. Goodbye."**

STEP 2: Enter Setting Mode

Start with the meter off. Then press and hold the **"SET"** button located on the front of the meter until the meter speaks. The meter is now in setting mode. See Figure 8.

Figure 8

Voice speaks: **"Your Prodigy® meter is on. You are now in setting mode."**

STEP 3: Setting Volume Level

Voice speaks: **"Volume Level 4."**

The meter displays **"VOL"** Press the **"Up"** or **"Down"** button to select the speaking volume. See Figure 9.

- Number 0 indicates that the speaking function is turned off.
- Numbers 1 to 7 indicate speaking volume from low to high, where "4" is displayed on the LCD during testing.


Figure 9

Then press the **"SET"** button and the year segment will appear and start flashing.

STEP 4: Set the Year

Voice speaks: **"The Year 2013."**

See Figure 10. Press the **"Up"** or **"Down"** button to advance the year. With the correct year on the display, press

the **"SET"** button and the date will appear on the display with the month segment flashing. See Figure 11.

STEP 5: Set the Month

Voice speaks: **"The Month November."**

Press the **"Up"** or **"Down"** button to advance the month. With the correct month on the display, press the **"SET"** button and the date segment will start flashing. See Figure 12.

STEP 6: Set the Date

Voice speaks: **"The Date 18."**

Press the **"Up"** or **"Down"** button to advance the date. With the correct date on the display, press the **"SET"** button and the time will appear on the display with the hour segment flashing. See Figure 13.

Figure 10

Figure 11

Figure 12

STEP 7: Set the Hour

Voice speaks: **"The Hour 10 PM."**

Press the **"Up"** or **"Down"** button to advance the hour. With the correct hour on the display, press the **"SET"** button and the minutes segment will start flashing. See Figure 14.

Figure 13


Figure 14

STEP 8: Set the Minutes

Voice speaks: **"The Minutes 08."**

Press the **"Up"** or **"Down"** button to advance the minutes. With the correct minutes on the display, press the **"SET"** button and the current unit of measurement will start flashing. See Figure 15

STEP 9: Select Unit of Measurement

Voice speaks: **"Blood glucose unit is mg/dL."**

Press and hold the **"Up"** or **"Down"** button until

the unit of measurement you are choosing appears on the display. Press the **"SET"** button and the memory segment will start flashing. Your meter can display test results in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L).

- The mg/dL unit is standard in the United States.
- The mmol/L unit is standard in Canada and Europe.

Important: Day averages are calculated from results obtained during the 7, 14, 21, 30, 60 and 90 days preceding the current date and time settings.

STEP 10: Delete Memory

When the **"dEL"** symbol and the **"M"** symbol appear on the display, you can choose to clear the memory. See Figure 16.

Voice speaks: **"Deleting record function: To skip, press 'Set' button. If you press the 'M' key for five (5) seconds, all records will be deleted."**


Figure 15

Figure 16

If you do not want to clear the memory, press the **"SET"** button to skip this step. If you want to clear ALL memory, press and hold the **"M"** button for five (5) seconds. **"- - -"** will appear to indicate all memory has been deleted. See Figure 17.


Figure 17

Setup is now completed. **"OFF"** is displayed before the meter shuts down. See Figure 18.

Figure 18

Voice speaks: **"Volume level one, the date November 18 2013. The time 10:08 pm. Milligrams per deciliter. Goodbye."**

Control Solution Testing

If your Prodigy® control solution did not come with a Prodigy® Control Solution Insert, the following information shall take the place of the insert.

Product Numbers:

990310 (Low)	53310 (Low)	53350 (High)
--------------	-------------	--------------

Intended Use:

Prodigy® control solution is a red liquid that contains glucose that will react with test strips and produce a test result. You can use either high or low Prodigy® control solution with your Prodigy® system.

Summary:

Prodigy® control solutions contain certain amounts of glucose that can be tested with all Prodigy® brand glucose meters. Test results should fall within the control ranges printed on the back of your Prodigy® test strip vial.

Caution: The control solution ranges are located on the back of the Prodigy® test strip vial. They are not recommended target ranges for your blood glucose.

Important Information:

- Use only Prodigy® control solutions with your Prodigy® meter.
- Check the expiration date on the bottle. Do not use if expired.
- Use within a period of 90 days from the date that you first open it. Record the discard date on the control solution bottle when you first open it to serve as a reminder to discard after 90 days.
- For in vitro diagnostic use.

Why Perform a Control Solution Test:

- To ensure that your meter and test strips are working properly together.
- To allow you to practice testing without using your own blood.

When to use:

- Once a week (to ensure that you continue to have accurate results).

- When you begin using a new vial of test strips.
- When test strips have been exposed to extreme environmental conditions.
- If you drop the meter.
- If you change the battery in the meter.

Storage and Handling:

- Store the control solution tightly closed at temperatures below 86°F (30°C).
- Do not refrigerate. Do not store in your car.

Composition:

	Glucose (% w/v)	Nonreactive Ingredients (% w/v)
Low Solution	0.05	99.95
High Solution	0.19	99.81

Performing a Control Solution Test

Start with the meter off.

STEP 1: Wash your Hands

Wash your hands with mild soap and water before performing any test. Be sure to dry them thoroughly.

STEP 2: Insert Test Strip

Insert a test strip with the contact bar end entering into the test strip slot first. Push the test strip as far as it will go without bending it. The meter turns on automatically.

STEP 3: Mark as a Control Solution Test

After the "d" symbol appears on the display, press the "M" button and a "M" symbol appears on the display. With the "M" symbol on the display, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the "M" button again and the "M" symbol will disappear.

Important: Be sure that you are in control solution mode so that the test result will not be stored in the meter memory.

STEP 4: Apply Control Solution

1. Check the expiration and discard dates on your control solution and test strip vials.
2. Shake the control solution bottle well, then remove cap.
3. Squeeze the bottle and discard the first drop then wipe the dispenser tip with clean tissue paper or cotton.
4. Squeeze the bottle again to get a second drop onto a clean, non-absorbent surface or on your clean fingertip. See Figure 19.
5. Bring the tip of the test strip to touch the drop of solution until the meter beeps.


Figure 19

Caution: To avoid contaminating the control solution with the content of the test strip, **DO NOT DIRECTLY APPLY CONTROL SOLUTION ONTO THE TEST STRIP.** See Figures 20 and 21.

Figure 20

Figure 21

STEP 5 : Check if the Test Result is in Range.

After the meter counts down from 6 to 1, the control solution test result appears on the display.

Compare the test result with the range printed on the test strip vial. Each vial of Prodigy® No Coding Test Strips may have a different control solution range. The result should fall within the printed range on the test strip vial.


Figure 22

See Figure 22.

Questionable Results:

If the test result falls outside the range printed on the test strip vial, check the "Trouble-Shooting Guide" located in the "System Troubleshooting" section and repeat the test.

Results that fall outside of the expected control ranges may indicate:

- An error in performing the test
- Control solution bottle not shaken well
- Expired or contaminated control solution
- Control solution that is too warm or too cold
- Test strip deterioration
- Meter malfunction
- Test strips expired

Caution: Do not use test strips or control solution that have exceeded the discard date, are expired or have been damaged. Your results may be inaccurate.

Contact Customer Care

If your test strips have been damaged or if you continue to get out-of-range results, it means that the system or the control solution may not be working properly. DO NOT use the system to test your blood glucose level. If you are unable to resolve the problem, contact Customer Care at **1.800.243.2636**.

Performing a Blood Glucose Test

Preparing a Lancing Device

STEP 1: Remove the cap of the lancing device by twisting it off. See Figure 23.


Figure 23

STEP 2: Insert a sterile lancet into the lancet holder and push down firmly until lancet is fully seated. Do not twist the lancet. See Figure 24.


Figure 24

STEP 3: Remove the protective cap from the lancet by twisting it and then save it for later use. See Figure 25.


Figure 25

STEP 4: Replace the cap onto the lancing device. Screw the cap until it is snug but not too tight.


Figure 26

STEP 5: Set the lancing level. The adjustable tip offers five (5) levels of skin penetration. To select the desired depth, twist the

adjustable tip in either direction until the number lines up with the arrow. Choose from the following depths: 1-2 (least penetration) for soft or thin skin, 3 (average penetration) for average skin, 4-5 (deepest penetration) for thick or calloused skin. See Figure 26.

Warning: To reduce the risk of infection:

- Never share a lancet or lancing device.
- Always use a new, sterile lancet.
- Lancets are for single use only, do not reuse.
- Avoid getting hand lotion, oils, dirt or debris in or on the lancets and the lancing device.
- These devices are intended to be used for patient self-monitoring and should not be used to collect blood from more than one person as this poses a risk of transmitting blood-borne pathogens such as Hepatitis B or HIV.

STEP 6: Cock the lancing device. Slide the cocking control back until it clicks. If it does not click, the lancing device may have been cocked when the lancet was inserted. See Figure 27.

Figure 27

*The lancing device is prepared and ready to lance your finger for a blood sample.

Getting a Blood Sample**STEP 1: Wash Your Hands and Puncture Site**

Use warm, soapy water. Rinse and dry your hands thoroughly.

STEP 2: Select and Lance a Puncture Site**• Fingertip**

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete. See Figure 28.


Figure 28

• Sites other Than Your Fingertip

Please refer to the "About Alternate Site Testing (AST)" Section. Please consult your healthcare professional before obtaining blood from sites other than your fingertip.

STEP 3: Obtain a Blood Sample

Do not smear the blood sample. See Figure 29. To obtain the most accurate results, wipe off the first drop of blood and gently milk the finger to obtain another drop of blood. See Figure 30.

Figure 29


Figure 30

STEP 4: Remove the Lancet

Remove the cap of the lancing device by twisting it off. Place the protective cap back on the exposed tip of the lancet by pushing the lancet into it. Eject the lancet by sliding the eject button forward. See Figure 31.

*Always use caution when removing the lancet. Discard the lancet according to your local safety regulations.


Figure 31

Warning: The first drop of blood usually contains tissue fluid and serum, which may affect the test result. It should be discarded.

Tips for Applying Blood Without Sight

Using the following four (4) tips may help you place a blood sample into the test strip successfully if you have little or no vision. The examples describe two (2) ways to bring the sample and the test strip together. Other techniques may work just as well. Keep experimenting until you find a technique that consistently works for you.

Tip One—Smaller has its Advantages

Because the Prodigy Voice® meter is so easy to move around, you can experiment with many different approaches to determine which works best for you. Get comfortable holding the meter in one hand while testing. The meter can be moved towards your finger as easily as your finger can be moved to the test strip.

Tip Two—Aim for the Tip

The tip of the test strip contains an opening to the absorbent channel in the center where the blood should be applied. The absorbent channel draws the blood into the test strip, much like a straw,

when the tip touches the blood sample. Blood should be applied only on the very tip. Blood applied to the bottom or sides of the test strip may cause inaccurate readings. It is important that the test strip does not bend at all when applying the sample.

Tip Three—Using a Sideways Approach

When testing, it may be helpful to rest or hold the meter on one side rather than flat or face up. In this position, the test strip sticks out of one end and the side with the buttons faces towards you.

The idea here is that the finger can be brought to the test strip from below until it touches the tip of the test strip. Then, the meter can be repositioned so the sample touches the absorbent opening at the end of the test strip. Notice how the meter can be positioned so that the test strip can be sticking out to the left or to the right. This allows for either hand to be used for testing.

Example 1

- Step 1—Place the meter on one side in front of you on a table with the test strip pointing to the left.

- Step 2—Gently squeeze a finger on the left hand and obtain a sample.
- Step 3—Bring your finger toward the tip of the test strip from the bottom. Place the puncture site on your finger to the tip of the test strip.
- Step 4—If you do not make contact with the blood and the meter does not begin the test, while keeping the finger in contact with the test strip, lift the rounded end of the meter from the table so the test strip is flush with the sample.
- Step 5—Move the test strip slightly, if necessary, until you hear the meter beep then say “now testing.” Hold the meter in place until you hear the test result. This procedure can be used with either hand by starting with the meter resting on the other side.

Example 2

- Step 1—Gently milk the finger to bring blood to the finger.
- Step 2—Puncture the finger.
- Step 3—Hold the meter in your hand. Place your index finger under the test strip to help

you guide the meter to where you have lanced.

- Step 4—Using overlapping dragging motions in an upward direction, repeat the motion around the puncture sight until you hear the meter beep then say, “now testing.”
- Step 5—Hold the test strip in place until you hear the meter say your test results.

Tip Four—Practice, Practice, Practice

As with all new skills, practicing is the key to feeling comfortable performing this new skill with confidence and ease. Don't give up; remember to practice, practice, practice. We suggest using Prodigy® control solution when practicing.

Practice until, with relative ease, you can get a test result that falls within the target range for the Prodigy® control solution three (3) tests in a row. It can be helpful to have someone with sight to watch you practice in order to provide feedback or to answer any questions you may have.

Applying Blood to the Test Strip

STEP 1: Insert the Test Strip

Insert a test strip with the contact bars end entering into the test strip slot first. Push the test strip as far as it will go without bending it. The meter turns on automatically. See Figure 32.

Figure 32

Voice speaks: **“Thank you for using Prodigy Voice®. Today’s date November 18, 2013. The time 10:08 am. Please apply blood to the test strip.”**

STEP 2: Apply Blood Sample

When the meter shows the “” symbol, apply blood to the opening of the absorbent channel of the test strip where it meets the narrow channel. Blood will be drawn into the test strip. See Figure 33.

Figure 33

Important: If you do not apply a blood sample within three (3) minutes, the meter will automatically turn off. You must remove the test strip and re-insert it again to turn on the meter and restart the test procedure.

The test strip confirmation window should be completely filled before the meter begins to count down. The meter will beep when the confirmation window is full. If the confirmation window is not full, “**L-b**” is displayed. You will have to use a new test strip.

Voice speaks: **“Now testing.”**

Caution: DO NOT:

- Smear or scrape the blood onto the test strip.
- Apply blood to the test strip when the test strip is not in the meter.
- Put blood or foreign objects into the test strip slot.

PRODIGY Voice® Blood Glucose Meter

If you have trouble filling the test strip, please contact Customer Care at **1.800.243.2636** for assistance.

STEP 3: View and Hear Your Result

After the meter counts down from 6 to 1, your blood glucose test result appears along with the unit of measurement.

Voice speaks: **"Your blood glucose 106 mg/dL."**

The blood glucose test result is automatically stored in the meter memory. Turn the meter off by removing the test strip. Discard the used test strip carefully to avoid contamination.

Figure 34

Voice speaks: **"Goodbye"**

Note: If at any time you need to hear the last message spoken, press the "Repeat" button and the last message will be repeated. See Figure 34.

PRODIGY Voice® Blood Glucose Meter

Caution: If you cannot test due to a problem with your testing supplies, contact Customer Care at **1.800.243.2636**. Failure to test could delay treatment and/or lead to serious medical conditions.

Using the Meter Memory

The Prodigy Voice® meter stores the 450 most recent blood glucose test results with date and time in its memory. It also provides you with 7, 14, 21, 30, 60 and 90-day averages of your blood glucose test results. You can review the individual or average test results.

OPTION 1: Recalling Individual Test Results

While the meter is turned off, press the “M” button to turn the meter on. After the prompts, press the “M” button again to enter memory mode and recall individual test results. Press the “Up” or “Down” buttons to advance to the next or previous test results. See Figure 35.


Figure 35

Note: Meter will not enter memory mode if test strip is inserted.

Voice speaks: **“Memory record for November 18, 2013, 10:08 PM. Your blood glucose 89 mg/dL.”**

Once you scroll through all the individual test results in the meter memory, “End” will display and the meter will turn off. See Figure 36.

Figure 36

Voice speaks: **“No memory record available. Goodbye.”**

Note: Once you are in the individual test recall mode, you cannot advance to the averaging simply by pressing the “M” button.

OPTION 2: Recalling Average Test Results

While the meter is off, press the “M” button and the meter will turn on.

Voice speaks: **“Your Prodigy® meter is on. Today’s date November 18, 2013; The time 10:08 PM. Please insert the test strip.”**


Figure 37

After the audible prompts, press and hold the “M” button for two (2) seconds and the meter will recall the 7-day average. See Figure 37.

Voice speaks: **"7-day average 106 mg/dL."**

Press the **"Up"** button to advance to the next average test result.

Note: The meter reports 7, 14, 21, 30, 60 and 90-day averages.

Exit the Memory Mode

To exit the memory mode during either option, press and hold the **"M"** button for approximately five (5) seconds and the meter will turn off.

Important: If you do not press any buttons for one (1) minute, the meter will display **"OFF"** and turn off automatically.

Specifications

Dimensions: 4.08 in.(L) x 2.31 in. (W) x 0.79 in. (H)
103.7 mm (L) x 58.8 mm (W) x 20 mm (H)

Weight: 3.10 oz with battery, 88g

Power Source: Two (2) 1.5 Volt AAA Alkaline Batteries

Battery Life: Over 1,000 tests

Display: LCD

Memory: 450 test results with date and time

External Output: USB Data Port

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Temperature warning

Operating Condition: 50°F - 104°F (10°C - 40°C),
below 85% R.H. (non-condensing)

Storage/Transportation Condition:

39.2°F - 104°F (4°C - 40°C), below 85% R.H.

Measurement Units: mg/dL or mmol/L

Measurement Range: 20 - 600 mg/dL (1.1 - 33.3 mmol/L)

The specifications may be changed without prior notice.

Viewing Results on a Personal Computer

Test results in your meter's memory can be transmitted to your personal computer. Please visit the Prodigy® website, www.prodigymeter.com for the latest instructions.

Caring for your Prodigy® Meter

Cleaning your System

Your Prodigy® meter should be cleaned whenever it is visibly dirty by wiping the outside of the meter using a cloth dampened with either mild detergent mixed with water or 70% isopropyl (rubbing) alcohol. Do not use bleach or other harsh abrasives to clean your meter.

Do not get any liquids, dirt, dust, blood or control solution inside the meter through the test strip port or the data port. Never spray cleaning solution on the meter or immerse it in any liquid.

Lancing Device

Wipe your Prodigy® lancing device with a soft cloth dampened with water and mild detergent or 70% isopropyl (rubbing) alcohol. Do not immerse the lancing device in any liquid. Do not use bleach.

Disinfecting your System

If your Prodigy® meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing

device should be decontaminated prior to use by the second person.

Disinfect your Prodigy® lancing device once a week. After cleaning the device, unscrew the cap and place it in 70% isopropyl (rubbing) alcohol for several minutes. Make sure the cap is completely dry before reattaching it to the device.

Please contact Customer Care if you have any questions at **1.800.243.2636**.

Storing Your System

Prodigy® meters do not require special maintenance.

- Store the meter, test strips and control solution in the carrying case after each use in a cool, dry place.
- Do not refrigerate.
- Do not store in the kitchen or bathroom.

*Your meter is a precision instrument. Please handle it with care.

Battery

The Prodigy Voice® meter comes with two (2) 1.5 V AAA alkaline batteries. Batteries are manufactured in a charged state and not designed for recharging. Recharging the batteries can cause battery leakage, or in some cases, high pressure rupture.

Replace both batteries at the same time.

Warning: Batteries can explode or leak and cause burns if installed backwards, disassembled, charged or exposed to water, fire or high temperatures.

Low Battery

The meter will alert you when the power is getting low by displaying two (2) different messages:

1. When the “” symbol appears alone on the display, the meter is functional and the results remain accurate, but you should change the batteries as soon as possible.

- When the “□” symbol appears together with the “E-b” symbol on the display, the batteries do not have enough power for a test. You must change the batteries before using the meter.

Battery Replacement

When replacing the batteries, make sure the meter is turned off.

STEP 1: With the meter off, press the buckle on the battery cover and lift up to open the cover.

STEP 2: Remove the old batteries and replace with new ones. Be sure to align the batteries properly. See Figure 38.

STEP 3: Close the battery cover. If the meter does not power on after you have replaced the batteries, check that the batteries are correctly installed.

STEP 4: Turn the meter on by pressing the “M” button to check the time and date. Replacing the


Figure 38

batteries does not affect the meter's memory (previous test results stored in the memory). However, the date, time and unit of measurement settings may need to be updated.

Caution: As with all small objects, the batteries should be kept away from small children. If the battery is swallowed, seek medical assistance immediately.

System Troubleshooting

Special Messages

Special messages and Error messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message.

In the event of a problem, refer to the information under “Action” in the “Error Messages” section. If you continue to have a problem, please refer to the “Troubleshooting Guide” section. If you follow the actions recommended but the problem is not resolved, please contact Customer Care at **1.800.243.2636** for assistance.

Message	What it means
	<p>“Lo” appears when your result is below the measurement limit, which is less than 20 mg/dL (1.1 mmol/L). “Lo” indicates hypoglycemia (low blood glucose). You should immediately consult your healthcare professional.</p>

Message	What it means
	<p>“HI” appears when your result is above the measurement limit, which is higher than 600 mg/dL (33.3 mmol/L). You should immediately consult your healthcare professional.</p>
Error Message	Description
	<p>What it means: Appears when the battery cannot provide enough power for a test. Action: Replace the battery immediately.</p>
	<p>What it means: Appears when inserting a used test strip. Action: Test with a new test strip. If the problem persists, please contact Customer Care at 1.800.243.2636.</p>

Error Message	Description
E-t	<p>What it means: Appears when the temperature is out of the system operating range (50°F-104° or 10°C-40°C).</p> <p>Action: Repeat the test after the meter and test strip are within the operating temperature range.</p>
L-b	<p>What it means: Insufficient blood sample.</p> <p>Action: Insert a new test strip.</p>

Only the most common error messages are listed. If your meter displays an error message that is not listed, please contact Prodigy® Customer Care at **1.800.243.2636**.

Troubleshooting Guide

- ❖ The meter does not display a message after inserting a test strip

Probable Cause	Actions
Battery exhausted.	Replace the batteries.
Battery incorrectly installed or absent.	Check that the batteries are correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the bar end entering into the test strip port first.
Defective meter.	Please contact Customer Care at 1.800.243.2636 for assistance.

- ❖ The test does not start after applying the sample.

Probable Cause	Actions
Insufficient blood sample.	Repeat the test using a new test strip with a larger blood sample.

PRODIGY Voice® Blood Glucose Meter

- ❖ The test does not start after applying the sample.

Probable Cause	Actions
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shut-off. (Two (2) minutes after user action).	Repeat the test with a new test strip. Apply sample only when the “  <p>Owner's Manual</p>

PRODIGY Voice® Blood Glucose Meter

- ❖ The control solution test is out of range.

Probable Cause	Actions
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Control solution bottle not shaken well.	Shake the control solution bottle vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date and the discard date of the control solution.
Control solution that is too warm or too cold.	Control solution should come to room temperature (less than 86°F/30°C) before testing.
Test strip deterioration.	Please repeat the test with a new test strip.

- ❖ The control solution test is out of range.

Probable Cause	Actions
Meter malfunction.	Please contact Customer Care at 1.800.243.2636 for assistance.

Information About your Prodigy® Meter

Comparing Meter and Laboratory Results

The test results you obtain from your meter may differ somewhat from your laboratory results due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way (see test strip package insert for typical accuracy and precision data as well as important information on limitations). To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before you go to the laboratory:

- Perform a control solution test to make sure that the meter is working properly.
- It is strongly recommended to fast for at least eight (8) hours before doing comparison tests.
- Take your meter with you to the laboratory.

While at the laboratory:

- Make sure that the samples for both tests (the meter test and the laboratory test) are

taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication or experienced stress.¹ In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test.²

Therefore, it is best to fast for eight (8) hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of bodily fluids (severe dehydration) may also cause a meter result to be different from a laboratory result.

References: ¹Surwit, R.S., and Feinglos, M.N.:

Diabetes Forecast (1988), April, 49-51.

²Sacks, D.B.: "Carbohydrates." Burtis, C.A, and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

Performance Characteristics

- **Accuracy:** ± 15 mg/dL when glucose < 75 mg/dL
 $\pm 20\%$ when glucose > 75 mg/dL
- **Precision:** This study shows the CV (correlation variation) is less than 5%.
- **The device has been certified to meet the following standards:**
98/79/EC, IEC 60601-1, IEC 61010-1, IEC 60601-1-2, IEC 61326, and ISO 15197

Expected Test Results

Please consult your physician to establish your target ranges. The American Diabetes Association (ADA) suggests the following targets for most nonpregnant adults with diabetes. More or less stringent goals may be appropriate for each individual.

Glycemic control Time of Day	ADA suggested reading
Preprandial plasma glucose (before a meal)	70-130 mg/dl
Postprandial plasma glucose (after a meal)	<180 mg/dl

Source: American Diabetes Association (2012. <http://www.diabetes.org/living-with-diabetes/treatment-and-care/blood-glucose-control/checking-your-blood-glucose.html>)

Your Prodigy® meter is designed to help you and your healthcare professional manage your diabetes. You must always rely on your healthcare professional to interpret your test results and to decide how to treat your diabetes.

Unexpected Test Results

Unexpected test results can occur. When this happens, please refer to these cautions.

Low Blood Glucose Test Results

False low results may occur if you are severely dehydrated. If you think you are severely dehydrated, contact your healthcare professional immediately.

If your test result is lower than 70 mg/dL or “Lo” appears on your meter screen, this could be a result of hypoglycemia. Although this result could be due to a testing error, it is best to treat first and retest later. Contact your healthcare professional immediately. See Figure 39.

Figure 39

High Blood Glucose Test Results

If your test result is higher than 180 mg/dL, this could be a result of hyperglycemia. If you are uncertain about this test result, retest. If the result continues to be higher than 180 mg/dL, contact your healthcare professional immediately.

If “HI” appears on your meter screen, this could be a result of severe hyperglycemia (results higher than 600 mg/dL). Retest your blood glucose. If “HI” appears again, contact your healthcare professional immediately. See Figure 40.

Figure 40

Repeated Unexpected Test Results

If you continue to get unexpected test results, check your blood glucose monitoring system by performing a control solution test. Refer to “Performing a Control Solution Test” section in this manual.

Unusual Red Blood Cell Count

A hematocrit range (percentage of red blood cells in the blood) that is extremely high (above 55%) or low (below 30%) can also cause false readings.

Accuracy of Meter Measurements

FAQ: Can two (2) blood glucose meters produce different readings that are from the same blood sample? Does this make one reading more accurate than the other?

Answer: Yes, different meters can produce different readings using the same blood sample. No, this does not necessarily mean that one result is more accurate than the other.

The reasons for this are as follows:

1. The FDA recognized standard for glucose meter accuracy is ISO 15197: In Vitro Diagnostic Test Systems Requirements for Blood Glucose Monitoring Systems for Self Testing in Managing Diabetes Mellitus. This standard sets accuracy requirements that manufacturers must meet before being approved by the FDA to market their glucose meters.

2. The accuracy requirements for the blood glucose monitoring system—meters and test strips—as stated in ISO 15197 are:

Ninety-five percent (95%) of the individual glucose results shall fall within ± 15 mg/dL (0.83 mmol/L) of the results of the manufacturer's measurement procedure at glucose concentrations < 75 mg/dL (< 4.2 mmol/L) and within $\pm 20\%$ at glucose concentrations ≥ 75 mg/dL (≥ 4.2 mmol/L).

How does this answer the FAQ

The accuracy of the system is determined by a clinical study using 100 blood samples at interval glucose concentrations ranging from < 50 mg/dL to over 400 mg/dL measured using both the manufacturer's meter and a clinical blood glucose analyzer for comparison (reference measurement).

Example

If the manufacturer's reference measurement is 70 mg/dL this means that the meter measurements are considered accurate if they fall within ± 15 mg/dL of 70 mg/dL (from 55 mg/dL to 85 mg/dL).

If meter A produces a measurement of 60 mg/dL and meter B measures 80 mg/dL, then both are in the acceptable range and meet the accuracy requirement.

Each manufacturer has met the System Accuracy requirement as part of its testing to be cleared for marketing. In the System Accuracy, testing a distribution of blood glucose readings will result at the various glucose concentration levels for one manufacturer's system. This is the same for other manufacturer's systems. Because of the range of acceptable readings for the accuracy requirement, identical readings from two different meters on the same blood sample may not happen.

Symbols Information

Symbol	Referent
	Do not re-use. Single use only.
	Consult Operating Instructions.
	Keep away from sunlight.
	Keep dry.
	Temperature limitation.
	Use by.
	Date of manufacture.
LOT	Batch code.
REF	Catalog number.
SN	Serial number.
CONTROL	Control.

Warranty Information

Limited Lifetime Warranty

Prodigy Diabetes Care, LLC extends a limited lifetime warranty to consumers who purchase a new Prodigy® brand meter. Under this limited lifetime warranty, your new meter is covered from the original date of purchase, as long as it has not been modified, altered, or misused. Under this limited lifetime warranty, Prodigy® will replace the meter free of charge if it is defective in material or workmanship.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE. PRODIGY DIABETES CARE, LLC WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING DIRECTLY OR INDIRECTLY FROM THE FAILURE OF THE PRODUCT TO PERFORM IN ACCORDANCE WITH SPECIFICATIONS.

Some states do not allow the exclusion or limitation of other express or implied warranties or incidental or consequential damages, so the above limitations or exclusions may not apply to you.

PRODIGY Voice® Blood Glucose Meter

PRODIGY®

Prodigy Diabetes Care, LLC 1.800.366.5901

PO Box 481928
Charlotte, NC 28269

For **Customer Care**
call: **1.800.243.2636**

www.prodigymeter.com

©2013 Prodigy Diabetes Care, LLC. Prodigy® and Prodigy Voice® are registered trademarks of Prodigy Diabetes Care, LLC. All other trademarks are property of their respective owners. All rights reserved. VOM653 Rev.2 12/13

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