
INSTRUCTIONS FOR TOUCH TALKING ANALOG WATCH

TOUCH TALKING WATCH MODEL (FULL FUNCTION)

- 1. Touch once to talk "Time"
- 2. Press and hold 3 seconds to talk "Date"

HOW TO SPEAK TIME AND DATE

- 1. In normal time, touch the watch face once to activate voice chip to speak time
- 2. In normal time, press and hold on the watch face for 3 seconds to speak date.



HOW TO ENTER MODE

1. Press the "Mode Key" for 3 seconds to enter mode. Then you will hear the circulate each time you press the button.

"Chime off" --> "Alarm Hour Setting" --> "Alarm minute setting" --> "Year setting" --> Month Setting --> "Date Setting" --> "12/24 Hour Setting" --> Time Hour Setting" --> Time Minute Setting --> "Beep Time Set" --> "Alarm is on" --> "Alarm is off"

2. Press the "Set Key" to do your setting.

HOW TO ENTER SLEEPING MODE

- 1. Press and hold the mode key for 6-8 seconds, In the meanwhile you will hear "chime on/off", please do not leave your finger until you hear "Do Do Do'.
- 2. After you hear "Do Do Do" please leave your finger for 1 second. Then press and hold the mode key again for 3 seconds to enter sleeping mode, you will hear a "beep beep" sound.
- 3. To wake up the watches, press and hold the mode key for 3 seconds again you will hear a "beep" sound.

CHECK ALARM STATUS

1. You can also check the alarm status by pressing the set key once. You will hear e.g. "Alarm is on 12 o'clock A.M."

** Please DO NOT touch the watch face when you are pressing the button to avoid duplicate function.